

Handicare 1100 Easy Hinge

A unique solution



Easy! Simply fold up the rail to clear the doorway

The lightweight Easy Hinge is effortlessly folded out of the way to clear the hallway while the stairlift is moved into a parked position. An ideal and cost-effective solution for hallways where access is required through doorways close to the staircase.

The folded rail also means there is no tripping hazard in a through corridor at the bottom of the stairs.

The Easy Hinge has a built-in gas spring, designed to make manual operation as easy as possible

Ideal for these situations:

- ✓ Doorway at the bottom of the staircase
- ✓ Open corridor or hallway
- ✓ Walking with rollator through the hallway
- ✓ Access with a wheelchair
- ✓ Lack of parking space for the stairlift



The Handicare 1100 can be parked anywhere on the rail. Even when it's parked on the hinged part the batteries will charge thanks to the continuous charge function throughout the rail.

Your local Handicare Business Partner



📍 Handicare Stairlifts B.V.
Newtonstraat 35
1704 SB Heerhugowaard
The Netherlands

✉️ stairlifts@handicare.com

🌐 www.handicare-stairlifts.com